

AMBU SPUR II RESUSCITATOR - infant



Product Code	34237
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

EAN/UPC	5707480032723
---------	---------------

Description

AMBU SPUR II INFANT - with face mask n°0

Ambu SPUR II provides users with exceptional tactile and visual feedback during resuscitation. The bag is highly responsive, with minimal mechanical resistance. Furthermore the characteristic design provides optimum stroke volume with perfect recoil.
Range includes infant, pediatric, and adult sizes.
All Ambu SPUR II resuscitators come in individual, resealable plastic carrying bags, complete with mask, oxygen tube and pressure limiting valve 4.0 kPa H2O.

The bags are color-coded for fast identification of infant, pediatric, and adult sizes.

Technical Specifications

Suitable for: infant < 10 kg
Stroke volume 1 hand 150 ml
Resuscitator volume: 220 ml
Size (length x diameter): 190x71 mm
Bag reservoir volume: 300 ml