

PVC SINGLE USE RESUSCITATOR KIT - adult



Product Code	34239
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279342390

Description

PVC RESUSCITATOR KIT - adult

Supplied with single use mask N° 5. Volume: 1,600 ml, oxygen reservoir 2,500 ml with valve, tongue holding forceps and oxygen tube in a nylon bag.

Disposable resuscitator bags with disposable mask, no return membrane valve and connection for PEEP valve.

PVC disposable, transparent mask shaped to allow a perfect adhesion to the patient's face. Latex free.