

# SILICONE RESUSCITATOR BAG with MASK N 5 - adult



Product Code	34260
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
--------	----

EAN/UPC	8023279342604
PARAF	965618729
GMDN	17591

## Description

### SILICONE AUTOCLAVABLE RESUSCITATORS - WITH HANDLE

- Type: Adult
- Mask N. 5
- Pressure relief valve H2O: 60 cm
- Volume ml: 1,600

Transparent resuscitators made of medical grade silicone and fully autoclavable at 121°C.  
Integrated handstrap for user comfort and complete control.

Made with no natural rubber or latex.

Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR.