

SILICONE RESUSCITATOR BAG with MASK N 3 - child



Product Code	34261
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
EAN/UPC	8023279342611
PARAF	965618756
GMDN	17591

Description

SILICONE AUTOCLAVABLE RESUSCITATORS - WITH HANDLE

- Type: Child
- Mask N. 3
- Pressure relief valve H2O: 40 cm
- Volume ml: 550

Transparent resuscitators made of medical grade silicone and fully autoclavable at 121°C.
Integrated handstrap for user comfort and complete control.

Made with no natural rubber or latex.

Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR.