



Product Code	34263
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
--------	----

RDM (NSIS)	2893262
CND	R03020201
EAN/UPC	8023279342635
GMDN	17591

Description

SILICONE AUTOCLAVABLE RESUSCITATORS KIT - WITH HANDLE

- Type: Adult
- Mask N. 5
- Pressure relief valve H₂O: 60 cm
- Volume ml: 1,600
- Reservoir bag: 2,500 ml
- 2 m oxygen tube
- Reservoir valve

Transparent resuscitators made of medical grade silicone and fully autoclavable at 121°C.
Integrated handstrap for user comfort and complete control.

Made with no natural rubber or latex.

Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR.