

# OXYGEN RESERVOIR BAG 2500 ml



Product Code	34275
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2625022
CND	R03020201
EAN/UPC	8023279342758
UMDNS	13360
PARAF	903140123
GTIN	08023279342758

## Description

To be used with GIMA codes: 34244 - 34245 - 34277 - 34200 (need valve 34279), 34230 - 34240, AMBU, LAERDAL and TECNOGAZ resuscitators.