



Product Code	34302
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2535677
CND	Z12021003
EAN/UPC	8023279343021
UMDNS	12293
PARAF	903147231
GTIN	08023279343021
GMDN	15076

Description

STANDARD RECHARGEABLE SETS (HANDLE + BLADES)
 • Adult set (3 blades Mc INTOSH 2-3-4)

- Great features of the unique GIMA blades and handles
- standard 2.5 V vacuum bulb
 - light ergonomic, satin finish handle to reduce glare
 - Stainless Steel contacts guarantee long working life
 - easy sterilization ETO or steam
 - full range of Mc Intosh & Miller blades, latex free
 - ISO 7376 fittings allow use of blades with existing handles
 - set is packaged in an anti-shock case and in cardboard box
 - 2.5 V handles work with both alkaline batteries or rechargeable batteries
 - good light transmission of 3,500 LUX with 2.5 V handles
 - autoclavable up to 134° for 5 minutes (approx. 2,000 times)
 - non magnetic blades