

MC-INTOSH BLADE N.1 - babies



Product Code	34318
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837798
CND	Z12021099
EAN/UPC	8023279343182
PARAF	903147357
GTIN	08023279343182
GMDN	46827

Description

- Blade N° 1 Mc INTOSH - babies
Overall Length: 95 mm
Distal Width: 10 mm
- Great features of the unique GIMA blades and handles
- standard 2.5 V vacuum bulb
 - light ergonomic, satin finish handle to reduce glare
 - Stainless Steel contacts guarantee long working life
 - easy sterilization ETO or steam
 - full range of Mc Intosh & Miller blades, latex free
 - ISO 7376 fittings allow use of blades with existing handles
 - set is packaged in an anti-shock case and in cardboard box
 - 2.5 V handles work with both alkaline batteries or rechargeable batteries
 - good light transmission of 3,500 LUX with 2.5 V handles
 - autoclavable up to 134° for 5 minutes (approx. 2,000 times)
 - non magnetic blades