



Product Code	34320
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837798
CND	Z12021099
EAN/UPC	8023279343205
PARAF	903147371
GTIN	08023279343205
GMDN	46827

## Description

Blade N° 3 Mc INTOSH - adults  
Overall Length: 134 mm  
Distal Width: 14 mm

- Great features of the unique GIMA blades and handles
- standard 2.5 V vacuum bulb
  - light ergonomic, satin finish handle to reduce glare
  - Stainless Steel contacts guarantee long working life
  - easy sterilization ETO or steam
  - full range of Mc Intosh & Miller blades, latex free
  - ISO 7376 fittings allow use of blades with existing handles
  - set is packaged in an anti-shock case and in cardboard box
  - 2.5 V handles work with both alkaline batteries or rechargeable batteries
  - good light transmission of 3,500 LUX with 2.5 V handles
  - autoclavable up to 134° for 5 minutes (approx. 2,000 times)
  - non magnetic blades