

MAXLITE F.O.LED ADULT SET 3 BLADES MC-INTOSH 2

-3-4 - 2.5 V - 20.000 lux



Product Code	34335
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2837787
CND	Z12021003
EAN/UPC	8023279343359
GTIN	08023279343359
GMDN	15076

Description

MAXLITE F.O. LED SET - 2.5 V

- Adult set (3 blades Mc Intosh 2-3-4)
- F.O. laryngoscopes provide cool bright illumination
- MaxliTE F.O. blades with integrated fibre
- LED F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades