

MAXLITE F.O.LED PEDIATRIC SET 3 BLADES MILLER 00-0-1 - 2.5 V - 20.000 lux



| | |
|---------------|----------------|
| Product Code | 34336 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 283778 |
| CND | Z12021003 |
| EAN/UPC | 8023279343366 |
| GMDN | 15076 |

Description

MAXLITE F.O. LED SET - 2.5 V

- Pediatric set (3 blades Miller 00-0-1)

- F.O. laryngoscopes provide cool bright illumination
- Maxlite F.O. blades with integrated fibre
- LED F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades