



Product Code	34340
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1831045
CND	Z1203020408
EAN/UPC	8023279343403
PARAF	903147460
GTIN	08023279343403
GMDN	45607

## Description

Small and lightweight wrist Oximeter available for home and clinical long-term and comfortable monitoring. It can be simply worn on the patient's wrist as a watch to start the monitoring and can be used for sports to measure SpO2 and pulse rate (PR). Not suitable for pediatric use.

- audible and visual alarm with adjustable alarm limits
- trend review function of SpO2 and PR
- selectable storage modes: 1 group/1s, 1 group/4s, and 1 group/8s
- up to 500 hours storage for SpO2 and PR in default mode (1 group/8s)
- clock setting and date display function
- integrated software to download and store all measurements on PC
- removable finger probe
- comfortable to be worn on wrist as a watch