



Product Code	34341
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
RDM (NSIS)	1698835
CND	Z1203020408
EAN/UPC	8023279343410
PARAF	979382963
GMDN	17148

Description

OXY 110 PULSE OXIMETER

- large 3.5 inch auto rotation TFT display for SpO2, PI, Plethysmogram, PR and temperature
- patented MoveOxy™ SpO2 technology with artifact removal and anti-motion performance
- various choice of SpO2 probes for adult, pediatric and neonate
- built-in lithium battery up to 18 hours continuous monitoring
- adjustable alarm limit, pitch tone level, key tone and brightness
- up to 580 hours SpO2 data storage
- upgradeable to infrared ear temperature

Supplied with cover, adult clip probe, charging base, USB cable.

Multilingual box and manual GB, FR, IT, ES

Technical Specifications

Technique: dual wavelength LED
 SpO2 range: 70 to 100% $\pm 2\%$
 Pulse rate range: 30 to 250 bpm $\pm 2\%$
 Perfusion index (PI) range: 0.2 to 20%
 Display: 3.5" TFT display
 Data storage: up to 580 hours
 Measuring mode: continuous or spot check
 Power supply: lithium battery
 Charging: charging stand or US B cable
 Dimension: 158x73x25 mm
 Weight: 230 g (including battery)
 Data upload: US B cable or wireless