

TRI-BALL RESPIRATORY EXERCISER



Product Code	33436
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	304238
CND	y030327
EAN/UPC	8436020761836
GMDN	31266

Description

TRI-BALL® - respiratory exerciser - box of 12

Tri-Ball® helps the patient to recover the normal respiration after a chest or abdominal surgery.
Easy to handle, light and break-resistant design.
Germ-proof mouthpiece holder.
It includes a particle filter.
Label identification to write patient name and room.
Single patient.

Label and user manual in GB, FR, DE, IT, ES, PT, NL, HU, PL, LT, CZ, SI, Arabic.

Made in Spain.