



Product Code	34361
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
EAN/UPC	4053755129157

## Description

### HEINE FIBRE OPTIC BLADES

- HEINE Code: F-000.22.763
- Blade type: XP MAC 3 - disposable (box of 25)
- Overall length: 134 mm
- Distal width: 14 mm

The "Classic+" range of F.O. blades with HEINE HiLite®

Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative HEINE LEDHQ light in the handle or the traditional XHL Xenon Halogen Bulb, the exclusive HEINE HiLite® Fibre Optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

#### Long lasting bulbs

HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! XHL 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. HEINE LEDHQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

#### Long lasting 3.5 V rechargeable batteries

HEINE 3.5 V systems (LEDHQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).