

HEINE CLASSIC+ LED LARYNGOSCOPE SET with 3.5V recharg.D6741 handle - 3 blades



Product Code	34367
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	1726137
CND	R9099
EAN/UPC	8023279343670

Description

HEINE CLASSIC+® FIBRE OPTIC LARYNGOSCOPE SETS with 3.5V rechargeable handle

- 3 blades set with Paed 1, Mac 2, Mac 3, LED HQ rechargeable handle
- HEINE Code: F-119.94.865

The "classic+" range of fibre optic blades with Heine Hilite®

Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative Heine led HQ light in the handle or the traditional XHI Xenon Halogen bulb, the exclusive Heine Hilite® fibre optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

5-years guarantee or 4,000 standard autoclave cycles of 134°C/5 min. with more than 1,000 lux on all blades with Heine Hilite® fibre bundle.

Long lasting bulbs Heine XHI® Xenon Halogen bulbs incorporate the latest Xenon Halogen technology to give you up to 40% more light! Heine XHI Xenon Halogen 3.5 V bulbs

last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. Heine led HQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

Long lasting 3.5 V rechargeable batteries Heine 3.5 V systems (led HQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery which avoids the need to replace non-rechargeable batteries at frequent intervals. Heine F.O. laryngoscope blades and handles are compatible with the iso 7376 (green standard).