



Product Code	34422
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2789843
CND	Z12021003
EAN/UPC	4053755127740
UMDNS	15-076

Description

HEINE FIBRE OPTIC BLADES

- HEINE Code: F-000.22.105
- Blade type: MAC 5
- Overall length: 176 mm
- Distal width: 15 mm

The “Classic+” range of F.O. blades with HEINE HiLite®
Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative HEINE LEDHQ light in the handle or the traditional XHL Xenon Halogen Bulb, the exclusive HEINE HiLite® Fibre Optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

5-years guarantee or 4,000 standard autoclave cycles of 134°C/5 min. with more than 1,000 lux on all blades with HEINE HiLite® fibre bundle.

Long lasting bulbs
HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! XHL 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. HEINE LEDHQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

Long lasting 3.5 V rechargeable batteries
HEINE 3.5 V systems (LEDHQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).