


REUSABLE LARYNGEAL MASK 2 - infant



|               |                |
|---------------|----------------|
| Product Code  | 34425          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |                |
|------------|----------------|
| UK-REP     | No             |
| CH-REP     | No             |
| RDM (NSIS) | 2526930        |
| CND        | R01020101      |
| EAN/UPC    | 8023279344257  |
| UMDNS      | 12448          |
| GTIN       | 08023279344257 |
| GMDN       | 45035          |

Description

SILICONE REUSABLE LARYNGEAL AIRWAY MASK



- 1. Optimal anatomical design ensure convenient and safe insertion. Universal 15 mm connector (ISO).
- 2. See-through fl exible airway tube, clear guideline for easy insertion.
- 3. Tiny protruding bumps minimize the vomiting substance fl owing from exophagus and into the lung.
- 4. Seamless cuff formation offers maximum comfort.
- 5. Strengthened oval tip prevents possible bending during insertion.
- 6. Thin walled pilot balloon facilitates the sensing of the infl ation tube.

Each laryngeal airway mask is individually serialized for identification and traceability purposes.

AUTOCLAVABLE UP TO 40 TIME.

Technical Specifications

- Size: 2
- Patient range: infant
- Weight range: 10-20 Kg
- Cuff volume (ml. of air): 10 ml