


REUSABLE LARYNGEAL MASK 2.5 - child



Product Code	34426
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2526930
CND	R01020101
EAN/UPC	8023279344264
UMDNS	12448
GTIN	08023279344264
GMDN	45035

Description

SILICONE REUSABLE LARYNGEAL AIRWAY MASK

- 1. Optimal anatomical design ensure convenient and safe insertion. Universal 15 mm connector (ISO).
- 2. See-through fl exible airway tube, clear guideline for easy insertion.
- 3. Tiny protruding bumps minimize the vomiting substance fl owing from exophagus and into the lung.
- 4. Seamless cuff formation offers maximum comfort.
- 5. Strengthened oval tip prevents possible bending during insertion.
- 6. Thin walled pilot balloon facilitates the sensing of the infl ation tube.

Each laryngeal airway mask is individually serialized for identification and traceability purposes.

AUTOCLAVABLE UP TO 40 TIME.

Technical Specifications

- Size: 2.5
- Patient range: child
- Weight range: 20-30 Kg
- Cuff volume (ml. of air): 14 ml