

GUEDEL AIRWAY 70 mm - boy - 1 - white



Product Code	34434
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2645421
CND	R010102
EAN/UPC	8023279344349
UMDNS	12732
PARAF	988194534
GTIN	08023279344349
GMDN	42424

Description

GUEDEL AIRWAYS STERILE - FLEXIBLE, TRANSPARENT

- Individually packed guedel airways
- Guedel airway 70 mm
- white
- boy - 1

Guedel airways in 9 sizes with smooth surface and round tip for improved patient comfort.
Made in polyethylene (PE). Latex free.