

GUEDEL AIRWAY 100 mm - adult medium - 4 - red



Product Code	34437
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2645421
CND	R010102
EAN/UPC	8023279344370
UMDNS	12732
GTIN	08023279344370
GMDN	32655

Description

GUEDEL AIRWAYS STERILE - FLEXIBLE, TRANSPARENT

- Individually packed guedel airways
- Guedel airway 100 mm
- red
- adult/medium - 4

Guedel airways in 9 sizes with smooth surface and round tip for improved patient comfort.
Made in polyethylene (PE). Latex free.