

HEINE RE-CHARGEABLE LI-ION L BATTERY X-007.99.383 - spare



Product Code	34448
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Yes
--------	-----

EAN/UPC	8023279344486
---------	---------------

Description

HEINE RECHARGEABLE LARYNGOSCOPE HANDLES AND BATTERIES

- Re-chargeable 3.5 V battery for codes 34372-3, 34402, 34367
- HEINE Code: X-007.99.383

The "classic+" range of fibre optic blades with Heine Hilite®

Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative Heine led HQ light in the handle or the traditional XHI Xenon Halogen bulb, the exclusive Heine Hilite® fibre optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

5-years guarantee or 4,000 standard autoclave cycles of 134°C/5 min. with more than 1,000 lux on all blades with Heine Hilite® fibre bundle.

Long lasting bulbs Heine XHI® Xenon Halogen bulbs incorporate the latest Xenon Halogen technology to give you up to 40% more light! Heine XHI Xenon Halogen 3.5 V bulbs

last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. Heine led HQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

Long lasting 3.5 V rechargeable batteries Heine 3.5 V systems (led HQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery which avoids the need to replace non-rechargeable batteries at frequent intervals. Heine F.O. laryngoscope blades and handles are compatible with the iso 7376 (green standard).