

"GIMA GREEN" PAEDIATRIC SET 2 BLADES MILLER 0-1



Product Code	34451
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2837804
CND	Z12021003
EAN/UPC	8023279344516
UMDNS	12293
PARAF	903148258
GMDN	15076

Description

GIMA GREEN F.O. XENON/HALOGEN SET - 2.5 V

- Pediatric set (2 blades Miller 0-1)
- F.O. laryngoscopes provide cool bright illumination
- Gima Green F.O. blades with detachable fibre
- Xenon/Halogen F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades