

"GIMA GREEN" PAEDIATRIC SET 2 BLADES MILLER 0-1



| | |
|---------------|----------------|
| Product Code | 34451 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2837804 |
| CND | Z12021003 |
| EAN/UPC | 8023279344516 |
| UMDNS | 12293 |
| PARAF | 903148258 |
| GTIN | 08023279344516 |
| GMDN | 15076 |

Description

GIMA GREEN F.O. XENON/HALOGEN SET - 2.5 V

- Paediatric set (2 blades Miller 0-1)
- F.O. laryngoscopes provide cool bright illumination
- Gima Green F.O. blades with detachable fibre
- Xenon/Halogen F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades