

"GIMA GREEN" PAEDIATRIC SET 3 BLADES MILLER 0-1, MC-INT.0



Product Code	34454
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2837804
CND	Z12021003
EAN/UPC	8023279344547
PARAF	926570704
GTIN	08023279344547
GMDN	15076

Description

GIMA GREEN F.O. XENON/HALOGEN SET - 2.5 V

- Pediatric set (3 blades Miller 0-1, Mc Intosh 0)

- F.O. laryngoscopes provide cool bright illumination
- Gima Green F.O. blades with detachable fibre
- Xenon/Halogen F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades