

USB GREEN FO.LARYNGOSCOPE HANDLE



Product Code	34468
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2219512
CND	Z12021003
EAN/UPC	8023279344684

Description

Universal USB handle for all green fibre optic blades.

This handle is a charging unit including lithiumion rechargeable battery and USB port for easy recharging without charging station.

Features:

- glass bead blasted, matt chromium plated surface for a high quality optic
- no surface structure, lateral grooves for hygienic use and secure grip in the hand
- metallic lamp case
- quick and easy removal of bulb
- autoclavable up to 134° (without batteries)
- handle in accordance with ISO 7376 (green standard)

Handle supplied in a rigid case with 3 adapters (US, GB, EU).

Made in Germany.