

# "GIMA GREEN" F.O. BLADE MC-INTOSH N. 5 - adult X-large



Product Code	34474
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837768
CND	Z12021003
EAN/UPC	8023279344745
PARAF	925813848
GTIN	08023279344745
GMDN	46827

## Description

"GIMA GREEN" FIBRE OPTIC BLADES - equipped with detachable fibre optic light guide

Mc INTOSH N° 5 blade - Adult x-large

### GIMA GREEN LARYNGOSCOPES

- maximum light transmission thanks to durable and high quality fibres
- specially designed blades for easy intubation
- matt finished to reduce glare
- autoclavable up to 134°C for 5 minutes. Blades autoclavable up to 2,000 times
- compatible to plasma sterilization
- all blades are compatible to BS EN ISO 7376 Green standard handles
- all blades are made from Austenitic Stainless steel thereby making them rust proof
- long life blades. GIMA offers 3 years Warranty
- a large selection of power sources to use: 'AA' and 'C' type batteries or 2.5 V and 3.5 V rechargeable batteries
- 2.5 V handles can be used with standard or rechargeable batteries
- handles are available with 2.5 V and 3.5 V Xenon/Halogen and LED bulbs
- all blades are non magnetic