

"GIMA GREEN" PAEDIATRIC HANDLE 2.5V



Product Code	34480
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2837774
CND	Z12021099
EAN/UPC	8023279344806
PARAF	903148322
GTIN	08023279344806
GMDN	46830

Description

Pediatric rechargeable metal handle (2xAA) - 2.5 V

F.O. XENON/HALOGEN BATTERY HANDLES - 2.5 V

- Closed battery compartment prevent liquid penetrating inside the battery
- Handles are supplied with lamp. Battery not included
- Solid metal handle of chromium plated metal

GIMA GREEN LARYNGOSCOPES

- maximum light transmission thanks to durable and high quality fibres
- specially designed blades for easy intubation
- matt finished to reduce glare
- autoclavable up to 134°C for 5 minutes. Blades autoclavable up to 2,000 times
- compatible to plasma sterilization
- all blades are compatible to BS EN ISO 7376 Green standard handles
- all blades are made from Austenitic Stainless steel thereby making them rust proof
- long life blades. GIMA offers 3 years Warranty
- a large selection of power sources to use: 'AA' and 'C' type batteries or 2.5 V and 3.5 V rechargeable batteries
- 2.5 V handles can be used with standard or rechargeable batteries
- handles are available with 2.5 V and 3.5 V Xenon/Halogen and LED bulbs
- all blades are non magnetic