

"GIMA GREEN" HALOGEN BULB 3.5V



Product Code 34491

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279344912

Description

Halogen Bulb 3.5 V

"GIMA GREEN" RECHARGEABLE HANDLES AND BATTERIES 3.5 V - OLD
up to lot 8 (ped handle) up to lot 7 (adult handle)

- Closed battery compartment prevent liquid penetrating inside the battery
- Handles are supplied with lamp. Battery not included
- Solid metal handle of chromium plated metal

GIMA GREEN LARYNGOSCOPES

- maximum light transmission thanks to durable and high quality fibres
- specially designed blades for easy intubation
- matt finished to reduce glare
- autoclavable up to 134°C for 5 minutes. Blades autoclavable up to 2,000 times
- compatible to plasma sterilization
- all blades are compatible to BS EN ISO 7376 Green standard handles
- all blades are made from Austenitic Stainless steel thereby making them rust proof
- long life blades. GIMA offers 3 years Warranty
- a large selection of power sources to use: 'AA' and 'C' type batteries or 2.5 V and 3.5 V rechargeable batteries
- 2.5 V handles can be used with standard or rechargeable batteries
- handles are available with 2.5 V and 3.5 V Xenon/Halogen and LED bulbs
- all blades are non magnetic