

"MAXLITE" FLEXI F.O. SET LED 2 BLADES MC-INTOSH 3, 4



| | |
|---------------|----------------|
| Product Code | 34498 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2837787 |
| CND | Z12021003 |
| EAN/UPC | 8023279344981 |
| GTIN | 08023279344981 |
| GMDN | 15076 |

Description

Deluxe set including one F.O. rechargeable handle, two F.O. Mc Intosh blades N° 3 and 4 and bulb.

The movable blade tip facilitates lifting of the epiglottis and enables a perfect visibility of the intubation process even in case of difficult anatomical conditions.

A better visibility and a lower lever force reduce the risk of injury of the larynx and the teeth.

- stainless steel construction, matt finish
- can be used with all Green Spec Handles
- compatible with ISO 7376 Green Standards
- channel provides improved viewing area
- 4.5 mm fibre bundle with more than 6,500 individual fibre strands
- high quality fibres ensuring excellent light transmission, 20,000 lux with 2.5V LED lamp
- 2 batteries type C or rechargeable batteries (not included)
- autoclavable at 134° C, approximately 2,000 times
- can sustain Plasma sterilization (Sterrad 100 NX)

Technical Specifications

- Fiber light carrier: integrated
- F.O. bundle Ø: 4.5 mm
- Individual fiber strands: > 6.500
- Light transmission with 2.5V: 6.000 lux
- Light transmission with 3.5V: 13.000 lux
- Autoclavable at 134 °C: 4000 times