

# OXYGEN CYLINDER 0.5 l with reducer - UNI - aluminium - empty



Product Code	34501
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279345018
---------	---------------

## Description

CYLINDERS FOR OXYGEN THERAPY (empty)

- UNI
- Litres: 0.5
- Flow: fixed 6 l
- Size: Ø 70x350mm
- Weight kg: 2.75
- Aluminium oxygen cylinder with pressure reducer to be assembled (not integrated)

Steel cylinders for therapeutic oxygen at 200 bar pressure, with valve and safety lock.  
Suitable for home and hospital therapeutic use, especially in case of desultory oxygen supply.

Made in Italy.