



Product Code	34502
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No
RDM (NSIS)	4524
EAN/UPC	8023279345025
PARAF	903140200
GMDN	35300

Description

CYLINDERS FOR OXYGEN THERAPY (empty)

- UNI
- Litres: 2
- Flow: 0/15 l
- Size: Ø100x400mm
- Weight kg: 4.5
- with soft bag + cover

Steel cylinders for therapeutic oxygen at 200 bar pressure, with valve and safety lock.
Suitable for home and hospital therapeutic use, especially in case of desultory oxygen supply.

Supplied in soft bag with cover. Oxygen cylinders are provided with a pressure reducer with two-stage release valve (0-15 l/min) and high pressure manometer to indicate oxygen quantity in the cylinder and face mask.

Made in Italy.