

OXYGEN CYLINDER 7 l with reducer - UNI - empty



Product Code	34507
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	4529
EAN/UPC	8023279345070
PARAF	903140275
GMDN	35300

Description

CYLINDERS FOR OXYGEN THERAPY (empty)

- UNI
- Litres: 7
- Flow: 0/15 l
- Size: Ø140x730mm
- Weight kg: 11
- with plastic bag + cover

Steel cylinders for therapeutic oxygen at 200 bar pressure, with valve and safety lock.
Suitable for home and hospital therapeutic use, especially in case of desultory oxygen supply.

Supplied in a sturdy plastic bag with cover. Oxygen cylinders are provided with a pressure reducer with two-stage release valve (0-15 l/min) and high pressure manometer to indicate oxygen quantity in the cylinder and face mask.

Made in Italy.