

INTEGRATED REDUCER with valve - UNI for 2/3 litres cylinders



Product Code	34510
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	83992
CND	Z120309
EAN/UPC	8023279345100
PARAF	903140337
GMDN	35300

Description

Pressure reducer with valve for 2/3 litres

- UNI

Pressure reducer with valve and index flowmeter adjustable 0-12 l/min.

It includes two stages, release valve and high pressure manometer.

Thread cylinder connection: E25/E17 conical refill thread connection UNI 21.7 x 1/14"

Apparatus to be connected to a cylinder without valve.