

OXYGEN CYLINDER 0.5 l with reducer - NF - aluminium - empty



Product Code 34521

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279345216

Description

CYLINDERS FOR OXYGEN THERAPY (empty)

- NF
- Litres: 0.5
- Flow: fixed 6 l
- Size: Ø 70x350mm
- Weight kg: 2.75
- Aluminium oxygen cylinder with pressure reducer to be assembled (not integrated)

Steel cylinders for therapeutic oxygen at 200 bar pressure, with valve and safety lock.
Suitable for home and hospital therapeutic use, especially in case of desultory oxygen supply.

Made in Italy.