

# INTEGRATED REDUCER with valve - NF for 2/3 litres cylinders



Product Code	34530
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	83992
CND	Z120309
EAN/UPC	8023279345308
PARAF	903140376
GMDN	35300

## Description

Pressure reducer with valve for 2/3 litres  
\*Products not in stock, available in 15 days

• NF

Pressure reducer with valve and index flowmeter adjustable 0-12 l/min.  
It includes two stages, release valve and high pressure manometer.  
Thread cylinder connection: E25/E17 conical refill thread connection NF Ø 22.91 x 1.514.  
Apparatus to be connected to a cylinder without valve.