

QV-600 PROFESSIONAL VEIN FINDER



Product Code	23454
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2420323
CND	V9099
EAN/UPC	8023279234541
UMDNS	14346
PARAF	950833982

Description

QV-600 PROFESSIONAL VEIN FINDER with display

Portable vein finding devices that can project blood vessels image on skin surface accurately and timely.
Ergonomic design with an advanced infrared imaging technology to provide real-time accurate vein digital image on patient's skin.
Help doctors and nurses to easily find veins of obese or chemotherapy patients, hairy or dark skin people, people with varicose veins, fine blood vessels and vascular obstruction.

Can adapt 12 different colours to different ages, body shapes, skin colours, weights and various operating environments.

Suitable for adults, children and infants.
Both allow to observe the vascular thickness, length, curve and branch information to determine the best puncture location to avoid pain.
Also available optional table stand or mobile trolley.

QV-600 additional features
-display with switchable languages (GB, FR, IT, ES, DE, PL, RU)
-possibility to photograph and store up to 6,000 images to track and analyze pathologies.

Technical Specifications

Light: harmless infrared light
Infrared wave: 760-940 nm
Display: 2.4"
Image: 3D pseudo-colour, AI algorithm
Adjustable brightness: 6 levels
Detection distance: 15-25 cm +/- 0.3 mm
Detection depth: 0-10 mm
Noise: ≈30 dB
Rechargeable Li-ion battery: 4800mAh
Power supply: 5V 2.0A, 100V-240V - 50/60Hz
Size - weight 24x6xh 6.5 cm - 450 g

Standard accessories

Aluminium carrying case
Calibration card
Charging cable
Adapter
User manual GB, FR, IT, ES