

QV-500 PROFESSIONAL VEIN FINDER



Product Code	23456
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2335065
CND	V9099
EAN/UPC	8023279234565
UMDNS	14346
PARAF	950833970

Description

QV-500 PROFESSIONAL VEIN FINDER

Portable vein finding devices that can project blood vessels image on skin surface accurately and timely.
Ergonomic design with an advanced infrared imaging technology to provide real-time accurate vein digital image on patient's skin.
Help doctors and nurses to easily find veins of obese or chemotherapy patients, hairy or dark skin people, people with varicose veins, fine blood vessels and vascular obstruction.

Can adapt 7 different colours to different ages, body shapes, skin colours, weights and various operating environments.

Suitable for adults, children and infants.
Both allow to observe the vascular thickness, length, curve and branch information to determine the best puncture location to avoid pain.
Also available optional table stand or mobile trolley.

QV-600 additional features
-display with switchable languages (GB, FR, IT, ES, DE, PL, RU)
-possibility to photograph and store up to 6,000 images to track and analyze pathologies.

Technical Specifications

Light: harmless infrared light
Infrared wave: 760-940 nm
Display: -
Image: -
Adjustable brightness: 5 levels
Detection distance: 15-25 cm +/- 0.3 mm
Detection depth: 0-8 mm
Noise: ≈20 dB
Rechargeable Li-ion battery: 4800mAh
Power supply: 5V 2.0A, 100V-240V - 50/60Hz
Size - weight 20x6.2xh 5.5 cm - 280 g

Standard accessories

Aluminium carrying case
Calibration card
Charging cable
Adapter

