

# PRESSURE REDUCER with flowmeter and humidif. - BRITISH STD.



Product Code	34571
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	84047
CND	Z120309
EAN/UPC	8023279345711
PARAF	903140477
GMDN	35300

## Description

Pressure reducer with release valve and H.P. manometer with adjustable ball flowmeter 0-15 l/min and autoclavable humidifier apparatus to be connected with cylinder valve.