



Product Code	34606
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2567510
CND	Z12159004
EAN/UPC	6933257923256
UMDNS	12873
PARAF	950833766

Description

SILENT OXYGEN CONCENTRATOR 5 I

Medical grade, silent, small oxygen concentrators for home use.
Provide high oxygen concentration through nasal cannula.
Work continuously.

Supplied with filter, humidifier bottle, nasal cannula and manual GB, IT.
Available on gimaitaly.com: ES, SE, GR, TR, Arabic.

Technical Specifications

Flow range: 0.5 ~ 5 l/min
Oxygen concentration: 95.5 - 87 %
Output pressure: 70 kPa
Input power: 350 VA
Noise level (front - whole unit): =45 dBA
Size - weight: 40x30xh 55 cm - 17.5 kg
Packing size: 48x40xh 63 cm
Operating voltage: 230 V - 50 Hz
Optimal oxygen concentration: after 15 min.
Classification: Type BF, Class II equipment