



Product Code	34647
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2371604
CND	M0304010103
EAN/UPC	8033140280614
PARAF	977814185

Description

ELASTIC COMPRESSION BANDAGES - individually packed, with clip.

PREVIFORTE - strong compression

- medical prescription bandage

- compression bandage for phlebitis treatment

- support bandage with strong, constant compression, ideal for varicose veins, for the prevention and cure of sport and work injuries.

Suitable for distortions, light fractures and sprains

Technical Specifications

- Composition:

- Cotton 92%

- Polyamide -

- Polyurethane 8%

- Extensibility:

- Longitudinal 170%

- Transversal -