

# COHESIVE COLOUR BANDAGE 4 m x 6 cm - mixed colours



Product Code 34668

Unit of sale box of 10 pcs

Minimum order 1

Type Medical device

Class I

UK-REP No

CH-REP No

RDM (NSIS) 2269546

CND M0303010201

EAN/UPC 8033140281789

PARAF 977814324

## Description

Cohesive bandages through microfine latex dots:

- For fixing dressings
- Ideal for supporting joint treating light fractures
- This bandage does not require clips because it remains well-attached to itself

## Technical Specifications

- Composition:
  - 42% cotton
  - 29% viscose
  - 28% polyamide
  - 1% elasthan and latex
- Width: 6 cm
- Length: stretched length 4 m
- Extensibility: approx. 90%