

COHESIVE COLOUR BANDAGE 4 m x 8 cm - red



|               |                |
|---------------|----------------|
| Product Code  | 34671          |
| Unit of sale  | box of 10 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2269546       |
| CND        | M0303010201   |
| EAN/UPC    | 8033140281703 |
| PARAF      | 977814351     |

Description

- Cohesive bandages through microfine latex dots:
- For fixing dressings
  - Ideal for supporting joint treating light fractures
  - This bandage does not require clips because it remains well-attached to itself

Technical Specifications

- Composition:
  - 42% cotton
  - 29% viscose
  - 28% polyamide
  - 1% elastan and latex
- Width: 8 cm
- Lenght: stretched lenght 4 m
- Extensibility: approx. 90%