

COHESIVE COLOUR BANDAGE 4 m x 8 cm - mixed colours



Product Code	34673
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2269546
CND	M0303010201
EAN/UPC	8033140281796
PARAF	977814375

Description

- Cohesive bandages through microfine latex dots:
- For fixing dressings
 - Ideal for supporting joint treating light fractures
 - This bandage does not require clips because it remains well-attached to itself

Technical Specifications

- Composition:
 - 42% cotton
 - 29% viscose
 - 28% polyamide
 - 1% elastan and latex
- Width: 8 cm
- Lenght: streched lenght 4 m
- Extensibility: approx. 90%