

ELASTIC TUBULAR NETTING 6 for head, shoulder, pelvis, leg



Product Code	34706
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2271842
CND	M03030103
EAN/UPC	8033971550092
PARAF	977814514
GMDN	13369

Description

Surgifix® elasticated tubular netting

- Material: Poly (polyester 38%, polyamide 35%, latex 27%)
- Range: 6
- Net for head, shoulder, pelvis, leg
- Length under tension: 25 m
- Width: 56 mm (Different width between Poly and Zerolatex is due to different elasticity of materials)

Surgifix® is a unique elastic tubular netting, light, soft and breathable.

It is the ideal solution for gauze or dressing retention.

Its unique structure makes it suitable for any part of the body, thus assuring perfect dressing adhesion under ideal compression grade and offering unrivalled patient's comfort.

Surgifix can be cut at the requested length.

Two ranges: Poly (polyester 38%, polyamide 35%, latex 27%) and Zerolatex (polyamide 86%, elastomer 14%).

Multilingual box (GB, FR, IT, ES, DE) and instructions (GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic).

Made in Italy.