

ELASTIC TUBULAR NETTING 9 for shoulder, body, back



Product Code	34709
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2271886
CND	M03030103
EAN/UPC	8033971550122
PARAF	977814540
GMDN	13369

Description

Surgifix® elasticated tubular netting

- Material: Poly (polyester 38%, polyamide 35%, latex 27%)
- Range: 9
- Net for shoulder, body, back
- Length under tension: 25 m
- Width: 85 mm (Different width between Poly and Zerolatex is due to different elasticity of materials)

Surgifix® is a unique elastic tubular netting, light, soft and breathable.
It is the ideal solution for gauze or dressing retention.
Its unique structure makes it suitable for any part of the body, thus assuring perfect dressing adhesion under ideal compression grade and offering unrivalled patient's comfort.
Surgifix can be cut at the requested length.
Two ranges: Poly (polyester 38%, polyamide 35%, latex 27%) and Zerolatex (polyamide 86%, elastomer 14%).

Multilingual box (GB, FR, IT, ES, DE) and instructions (GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic).

Made in Italy.