

# COHESIVE NON WOVEN ELASTIC BANDAGE 4.5 m x 10 cm - blue



Product Code	34738
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2269553
CND	M0303010201
EAN/UPC	8033140283004

## Description

### PREVITNT NON WOVEN ELASTIC BANDAGES

For fixing dressings. Ideal for supporting joints, treating light fractures.  
This bandage does not require clips because it remains well-attached to itself.  
Particularly suitable for use on joints, for sports and veterinary use.

- Composition:
  - 85% polypropylene
  - 15% spandex and latex.
- Extensibility: approx 100%.