

COHESIVE NON WOVEN ELASTIC BANDAGE 4.5 m x 10 cm - 5 colours



Product Code	34740
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2269553
CND	M0303010201
EAN/UPC	8033140281918

Description

PREVITNT NON WOVEN ELASTIC BANDAGES
For fixing dressings. Ideal for supporting joints, treating light fractures.
This bandage does not require clips because it remains well-attached to itself.
Particularly suitable for use on joints, for sports and veterinary use.

Mix 2 bandages x 5 colours (red, green, blue, yellow, skin).

- Composition:
 - 85% polypropylene
 - 15% spandex and latex.
- Extensibility: approx 100%.