



| | |
|---------------|----------------|
| Product Code | 34745 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2450789 |
| CND | M03040299 |
| EAN/UPC | 8023279347456 |
| UMDNS | 10-274 |
| PARAF | 973923016 |
| GTIN | 08023279347456 |
| GMDN | 10284 |

Description

KINESIOLOGY STANDARD TAPES - thickness 0.48 mm

- Kinesiology tape 5 m x 5 cm - black

High quality Kinesiology tape conforms to all areas of the body. Provides firm support and compression. Ideal for treatment of muscular symptoms, stimulus to blood circulation and lymph and fascia diversion. Reduce pain improving athletic performance. Easy to tear and unwinds to the core.

Individually packed in a plastic bag and cardboard box.

- soft hypoallergenic cotton cloth, elasticity up to 160%
- 70% water proof (level 3) and sweat proof
- supersticky: up to 7 days

Multilingual box and instructions (GB, FR, IT, ES)

Made of 97% cotton, 3% spandex