



Product Code	34752
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2450789
CND	M03040299
EAN/UPC	8023279347524
UMDNS	10-274
PARAF	973923079
GTIN	08023279347524
GMDN	10284

## Description

KINESIOLOGY STANDARD TAPES - thickness 0.48 mm

• Kinesiology tapes 5 m x 5 cm - mix colours

Mixed colours box includes 6 rolls in different colours: ceil blue, red, pink, green, black, beige

High quality Kinesiology tape conforms to all areas of the body. Provides firm support and compression. Ideal for treatment of muscular symptoms, stimulus to blood circulation and lymph and fascia diversion. Reduce pain improving athletic performance. Easy to tear and unwinds to the core.

Individually packed in a plastic bag and cardboard box.

- soft hypoallergenic cotton cloth, elasticity up to 160%
- 70% water proof (level 3) and sweat proof
- supersticky: up to 7 days

Multilingual box and instructions (GB, FR, IT, ES)

Made of 97% cotton, 3% spandex